

pH ELECTRODES

All electrodes are supplied with a 1m lead and BNC connector. Other plug types (for older meters) and cable lengths are available on request.

COMBINATION ELECTRODES

	Type	Description	Temp. range, °C
C51-8010	Combination	Standard glass, general purpose	0 to 100
C51-8060	Combination	Gel-filled, plastic body, general purpose	0 to 80
C51-8230	Spear combination	Glass body, for soil measurements	0 to 100
C51-8260	Flat head combination	Glass body, for surface measurements	0 to 100
C51-8310	Combination	Glass body, long reach, thin stem	0 to 100
C51-8360	Combination	Glass body, semi-micro	0 to 100

pH GLASS ELECTRODES

For use with separate reference electrode

	Type	Description	Temp. range, °C
C51-8520	Glass (indicator)	For use with calomel reference electrode	0 to 60
C51-8560	Glass (indicator)	For use with Ag/AgCl reference electrode	0 to 100

REFERENCE ELECTRODES

	Type	Description	Temp. range, °C
C51-8710	Reference	Sealed cap, calomel 4mm connector	0 to 60
C51-8760	Reference	Sealed cap, double junction Ag/AgCl, 4mm connector	0 to 100

METAL (REDOX) ELECTRODE

	Type	Description
C51-8910	Metal combination	Platinum for ORP measurement

