

BUFFERS

Buffer Tablets

Each tablet makes 100ml buffer solution. Pack of 50 tablets.

C51-9610 Buffer pH 4.0

C51-9630 Buffer pH 7.0

C51-9650 Buffer pH 9.2

ION SELECTIVE ELECTRODES

For the detection of the activity of specific ions in aqueous solution. Combination electrodes therefore do not require reference electrodes. Supplied with BNC connector. Must be used with pH/ion meters having a concentration scale or mV scale reading to 0.2mV.

	Type	Concentration range
C51-9800	Ammonia	1 to 10^{-7} M
C51-9850	Chloride	1 to 10^{-5} M
C51-9880	Silver/sulphide	1 to 10^{-7} M (both ions)
C51-9920	Sodium	1 to 10^{-6} M

pH Indicator Papers - see Paper section



C51-9800 et seq