

PIPETTES**PIPETTES, Class B**

One-mark soda-lime glass pipettes to BS1583, Class B, colour coded.

	Capacity, ml	Colour code
C52-1000	1	Blue
C52-1020	2	Orange
C52-1040	3	Black
C52-1060	4	2 Red
C52-1080	5	White
C52-1100	10	Red
C52-1120	15	Green
C52-1140	20	Yellow
C52-1160	25	Blue
C52-1180	50	Red
C52-1200	100	Yellow

PIPETTES, Class A

One-mark soda-lime glass pipettes to BS1583, Class A, colour coded.

	Capacity, ml	Colour code
C52-1400	1	Blue
C52-1420	2	Orange
C52-1440	3	Black
C52-1460	4	2 Red
C52-1480	5	White
C52-1500	10	Red
C52-1520	15	Green
C52-1540	20	Yellow
C52-1560	25	Blue
C52-1580	50	Red
C52-1600	100	Yellow

PIPETTES, Plastic

One-mark pipettes in polypropylene.

	Capacity, ml
C52-3000	1
C52-3020	2
C52-3040	5
C52-3060	10
C52-3080	25
C52-3100	50



C52-1000 et seq



C52-1400 et seq



C52-3000 et seq

GRADUATED PIPETTES

GRADUATED PIPETTES, Type 1, Class B

Class B pipette in soda-lime glass with permanent amber markings. Zero at the top.

	Capacity x graduation, ml	Colour code
C52-3300	1 x 0.01	Yellow
C52-3320	2 x 0.02	Black
C52-3340	5 x 0.05	Red
C52-3360	10 x 0.10	Orange
C52-3380	25 x 0.20	Green

GRADUATED PIPETTES, Type 2, Class B

Class B pipette in soda-lime glass with permanent amber markings. Zero at the jet; not to be blown out.

	Capacity x graduation, ml	Colour code
C52-3500	1 x 0.01	Yellow
C52-3520	2 x 0.02	Black
C52-3540	5 x 0.05	Red
C52-3560	10 x 0.10	Orange
C52-3580	25 x 0.20	Green

GRADUATED PIPETTES, Type 3, Class B

Class B pipette in soda-lime glass with permanent amber markings. Zero at the top; not to be blown out.

	Capacity x graduation, ml	Colour code
C52-3800	1 x 0.01	Yellow
C52-3820	2 x 0.02	Black
C52-3840	5 x 0.05	Red
C52-3860	10 x 0.10	Orange
C52-3880	25 x 0.20	Green

GRADUATED PIPETTES, Type 1, Class A

Class A pipette in soda-lime glass with permanent amber markings. Zero at the top.

	Capacity x graduation, ml	Colour code
C52-4400	1 x 0.01	Yellow
C52-4420	2 x 0.02	Black
C52-4440	5 x 0.05	Red
C52-4460	10 x 0.10	Orange
C52-4480	25 x 0.20	Green

GRADUATED PIPETTES, Type 2, Class A

Class A pipette in soda-lime glass with permanent amber markings. Zero at the jet.

	Capacity x graduation, ml	Colour code
C52-4800	1 x 0.01	Yellow
C52-4820	2 x 0.02	Black
C52-4840	5 x 0.05	Red
C52-4860	10 x 0.10	Orange
C52-4880	25 x 0.20	Green



C52-3300 et seq



C52-3500 et seq



C52-3800 et seq

GRADUATED PIPETTES, PLASTIC

Graduated pipette, Class B, in polypropylene.

	Capacity x graduation, ml
C52-6900	1 x 0.1
C52-6920	2 x 0.1
C52-6940	5 x 0.1
C52-6960	10 x 0.1

SEROLOGICAL PIPETTES

SEROLOGICAL PIPETTES, DISPOSABLE, VOLAC

Disposable glass serological pipettes, sterile, individually wrapped.

	Capacity x graduation, ml	Pack of 250
C52-7400	1 x 0.01	
C52-7420	1 x 0.1	
C52-7440	2 x 0.02	
C52-7460	5 x 0.05	


C52-7480	10 x 0.1	Pack of 125
-----------------	----------	-------------

DISPOSABLE PIPETTES

DISPOSABLE PIPETTES, PLASTIC

Sterile, plastic disposable pipettes

	Capacity	Type	Sleeve qty.	Pack of 500
C52-7700	1	Taper jet	10	
C52-7720	1	Taper jet	1	
C52-7740	1	Open end	10	
C52-7760	2	Taper jet	25	
C52-7780	2	Taper jet	1	
				Pack of 200
C52-7800	5	Taper jet	10	
C52-7820	5	Taper jet	1	
C52-7840	5	Open end	10	
				Pack of 500
C52-7940	10	Open end	10	
C52-7960	10	Moulded jet, moulded adapter	25	
C52-7980	10	Moulded jet, moulded adapter	1	



C52-6900 et seq



C52-7400 et seq



C52-7700 et seq



C52-7960

PASTEUR PIPETTES

PASTEUR PIPETTES, GLASS

	Type	Sleeve qty.	Pack of 1000
C52-8200	Short form	250	
C52-8220	Long form	250	
C52-8240	Short form, plugged	250	
C52-8260	Long form, plugged	250	

PASTEUR PIPETTES, PLASTIC

Length 152mm, bulb can draw up 3.5ml with one squeeze.
Capacity indicates extent of graduations on stem, at 0.5ml intervals.
Sterile pipettes are supplied in inner sleeves of 20.

	Capacity	Sterility	
C52-8500	1ml	Non-sterile	Pack of 500
C52-8520	3ml	Non-sterile	Pack of 500
C52-8550	1ml	Sterile	Pack of 1000
C52-8570	3ml	Sterile	Pack of 1000

DROPPING PIPETTES

DROPPING PIPETTES, GLASS

Glass dropping pipette with teat, length approximately 75mm.

C52-9900 Dropping Pipette, with teat. Pack of 10.

DROPPING PIPETTES, PLASTIC

Polyethylene dropping pipette with integral bellows, pack of 10.

	Volume, ml	Stem length, mm
C53-0100	2.5	80
C53-0120	4	120

PIPETTE ACCESSORIES

BULBS, PEAR SHAPED

PVC, for use with pipettes etc. Pack of 10

	Capacity, ml	Neck i.d., mm	Bulb dia., mm
C19-5800	22	5	39
C19-5850	42	5	50
C19-5900	60	5	58
C19-5950	120	11	68
C19-6000	240	16	80

BULBS (TEATS)

Rubber for use with dropper pipettes. Pack of 20

	Capacity, ml	Neck id., mm
C19-6200	1.5	4
C19-6250	2.0	5
C19-6300	2.75	6
C19-6350	5.0	6



C52-8200 et seq




C52-9900 et seq



C53-0100



C19-5800 et seq



C19-6200 et seq