

PIPETTORS

PIPETTORS, FIXED VOLUME, VOLAC

High-precision fixed volume micropipettors which fit comfortably in the hand. Each model is colour coded and marked with capacity and serial number. With tip ejector button. Lower body can be removed for autoclaving.

	Volume, μl	Colour code	Accuracy within $\pm\%$
C53-2010	5	Purple	1.0
C53-2030	10	Grey	1.0
C53-2070	20	Lime	0.5
C53-2090	25	Blue	0.5
C53-2150	50	Yellow	0.5
C53-2270	100	aqua	0.4
C53-2330	200	Orange	0.3
C53-2350	250	Grey	0.3
C53-2410	500	Blue	0.3
C53-2490	1000	Purple	0.3
C53-2510	5000	Lime	0.3

PIPETTORS, VARIABLE VOLUME, DIGITAL, VOLAC

High precision, variable volume, digital micropipettors with new slim design. Each model is colour coded and marked with the volume range and serial number. With tip ejector button and instruction leaflet. Lower body can be removed for autoclaving.

	Range, μl	Colour code	Tested at, μl	Accuracy within $\pm\%$
C53-2710	5 to 50	Yellow	5	2.0
			25	1.0
			50	0.5
C53-2760	10 to 100	Aqua	10	1.5
			50	0.6
			100	0.4
C53-2810	20 to 200	Orange	20	1.0
			100	0.2
			200	0.3
C53-2860	50 to 500	Blue	50	1.5
			250	0.75
			500	0.4
C53-2910	100 to 1000	Purple	100	0.8
			500	0.6
			1000	0.4
C53-2960	500 to 5000	Lime	500	0.8
			2500	0.4
			5000	0.4

PIPETTOR RACK

Plastic rack to hold up to 5 VOLAC single channel pipettors.

C53-2980 Pipettor Rack



C53-2010 et seq



C53-2710 et seq

MULTICHANNEL PIPETTORS, VARIABLE VOLUME, VOLAC

8 channel, variable volume, air displacement pipettor. Light weight design, easy to maintain and recalibrate, with tip ejector. All models take a standard 200µl yellow tip, but model **C53-3040** will also take a 300µl white tip for the larger volumes.

	Range, µl	Tested at, µl	Accuracy within ±%
C53-3000	0.1 to 10	1	0.4
C53-3020	5 to 50	5	0.15
C53-3040	40 to 300	40	0.6

DISPOSABLE PIPETTE TIPS

PIPETTE TIPS, STANDARD

Polypropylene pipette tips, designed to fit Volac pipettors and many other leading brands. Autoclavable. Bulk packed.

	Capacity, µl	Tip colour	Pack size
C53-3090	0 - 20	White	1000
C53-3100	1 - 200	Yellow	1000
C53-3150	1 - 300	White	1000
C53-3200	100 - 1000	Blue	1000
C53-3250	500 - 5000	White	500

PIPETTE TIPS, RACKED

As **C53-3100** and **C53-3200** but packed vertically in open rack to allow one-handed fitting to the pipettor.

	Capacity, µl	Tip colour	Pack size
C53-3300	1 - 200	Yellow	1000 (5 x 200)
C53-3320	100 - 1000	Blue	1000 (10 x 100)

PIPETTE TIPS, RACKED AND BOXED

As **C53-3100** and **C53-3200** but packed vertically in a rack within a box, for storage or autoclaving.

	Capacity, µl	Tip colour	Pack size
C53-3360	1 - 200	Yellow	96
C53-3380	100 - 1000	Blue	100

PIPETTE TIPS, RACKED AND BOXED, STERILE

As **C53-3100** and **C53-3200** but packed vertically in a rack within a box, pre-sterilised.

	Capacity, µl	Tip colour	Pack size
C53-3410	1 - 200	Yellow	96
C53-3460	100 - 1000	Blue	100



C53-3000 et seq



C53-3360 et seq

STOP PRESS !!

NEW PRODUCT

Economy Range of Variable Volume Pipettors

A range of simple to use automatic pipettors, covering 0.1 - 10ml with an accuracy of $\pm 1.75\%$.

This is an ideal teaching pipettor as it replaces the pipette filler and glass pipette, offering much greater safety in the laboratory along with hands-on experience of micropipettor use.

Fully autoclavable, with removable tip ejector, easy volume selection by micrometer style thumb press.

C53-2600 Eco Pipettor, 0.1 - 1.0ml

C53-2610 Eco Pipettor, 0.5 - 5.0ml

C53-2620 Eco Pipettor, 1.0 - 10.0ml

Accessories

C53-2650 Tip, 1000 μ l, blue, pack of 50

C53-2660 Tip, 5000 μ l, white, pack of 25

C53-2670 Tip, 10000 μ l, blue, pack of 25

