

MULTICHANNEL PIPETTORS, VARIABLE VOLUME, VOLAC

8 channel, variable volume, air displacement pipettor. Light weight design, easy to maintain and recalibrate, with tip ejector. All models take a standard 200µl yellow tip, but model **C53-3040** will also take a 300µl white tip for the larger volumes.

	Range, µl	Tested at, µl	Accuracy within ±%
C53-3000	0.1 to 10	1	0.4
C53-3020	5 to 50	5	0.15
C53-3040	40 to 300	40	0.6

DISPOSABLE PIPETTE TIPS

PIPETTE TIPS, STANDARD

Polypropylene pipette tips, designed to fit Volac pipettors and many other leading brands. Autoclavable. Bulk packed.

	Capacity, µl	Tip colour	Pack size
C53-3090	0 - 20	White	1000
C53-3100	1 - 200	Yellow	1000
C53-3150	1 - 300	White	1000
C53-3200	100 - 1000	Blue	1000
C53-3250	500 - 5000	White	500

PIPETTE TIPS, RACKED

As **C53-3100** and **C53-3200** but packed vertically in open rack to allow one-handed fitting to the pipettor.

	Capacity, µl	Tip colour	Pack size
C53-3300	1 - 200	Yellow	1000 (5 x 200)
C53-3320	100 - 1000	Blue	1000 (10 x 100)

PIPETTE TIPS, RACKED AND BOXED

As **C53-3100** and **C53-3200** but packed vertically in a rack within a box, for storage or autoclaving.

	Capacity, µl	Tip colour	Pack size
C53-3360	1 - 200	Yellow	96
C53-3380	100 - 1000	Blue	100

PIPETTE TIPS, RACKED AND BOXED, STERILE

As **C53-3100** and **C53-3200** but packed vertically in a rack within a box, pre-sterilised.

	Capacity, µl	Tip colour	Pack size
C53-3410	1 - 200	Yellow	96
C53-3460	100 - 1000	Blue	100



C53-3000 et seq



C53-3360 et seq