

FACESHIELD, ACETATE

Faceshield, comprising browguard with adjustable headband and 200mm acetate visor giving full face protection from chemical splash and impact. Conforms to BS EN 166:2002.

C59-4800 Faceshield, Acetate

FACESHIELD, POLYCARBONATE

Faceshield, comprising browguard with adjustable headband and 200mm polycarbonate visor giving full face protection from mechanical impact. Conforms to BS EN 166:2002.

C59-4850 Faceshield, Polycarbonate

EYE WASH BOTTLE

White polyethylene bottle, capacity 500ml, contoured eye-cup on a curved stem, with fine washing spray, overflow tube and dust cap. Bottle supplied empty.

C59-5000 Eye Wash Bottle

EYE WASH POST

Eye wash bottle, capacity 500ml, and small eye dropper bottle on wall mounting polystyrene support, featuring printed instructions for use. Bottles are supplied empty.

C59-5100 Eye Wash Post, Economy

EYE WASH STATION

Wall mounting eye wash station with 500ml eye wash bottle, green refill bottle, and eye dropper bottle. Held in place on white polystyrene back board by large clips. Back board also features instructions for use and a mirror. Overall dimensions 400 x 400mm. Bottles are supplied empty.

C59-5150 Eye Wash Station

EYE WASH SOLUTION

Suitable for all above eye wash bottles, 0.9% saline solution, sterile.

C59-5250 Eye Wash Solution, pack of 2 x 500ml

EMERGENCY EYE WASH SPRAY UNIT

Manually operated eye wash with 1.5m flexible tube. Provides an aerated spray at low pressure for washing the face and/or eyes or body of the injured person. Intended for use in laboratories in industry, schools or universities. Requires a mains or emergency feed, preferably drinking water, minimum pressure 2 bar, 3/4" connection. Eye wash should be positioned in an easily accessible place, no more than 10 seconds away from any potential hazard.

C59-5300 Emergency Eye Wash Spray Unit



C59-4800 & C59-4850



C59-5000



C59-5100



C59-5150



C59-5300