

## FILTER TUBES, PYREX

For small scale filtration applications. Side arm is 6mm diameter.  
Pack of 10.

**C65-6400** Filter Tube, Pyrex, 125 x 16mm, length x tube o.d.

**C65-6420** Filter Tube, Pyrex, 150 x 18mm, length x tube o.d.

**C65-6440** Filter Tube, Pyrex, 150 x 24mm, length x tube o.d.

## TEST TUBES, PYREX, STOPPERED

Borosilicate glass tubes with ground socket and polyethylene stopper.  
Pack of 10.

size	Size, mm	Capacity	Stopper
<b>C65-6600</b>	100 x 13	5ml	10/19
<b>C65-6620</b>	125 x 18	10ml	14/23
<b>C65-6640</b>	150 x 22	25ml	19/26
<b>C65-6660</b>	200 x 26	50ml	24/29

## TEST TUBES, PYREX, GRADUATED, STOPPERED

Borosilicate glass, graduated tubes with ground socket and polyethylene stopper.

	Size, mm	Capacity x subdivision	Stopper size	Pack size
<b>C65-6900</b>	100 x 14	5ml x 0.1ml	10/19	10
<b>C65-6920</b>	125 x 18	10ml x 0.2ml	14/23	10
<b>C65-6940</b>	150 x 22	25ml x 0.5ml	19/26	5
<b>C65-6960</b>	200 x 26	50ml x 1.0ml	24/29	2

## DURHAM TUBES

For fermentation tests. To BS EN ISO 4142:2002. Pack of 250

**C65-7200** Durham Tube, Soda Glass, 25 x 6.5mm

**C65-7220** Durham Tube, Soda Glass, 30 x 6.5mm

**C65-7240** Durham Tube, Soda Glass, 35 x 8mm

**C65-7260** Durham Tube, Soda Glass, 50 x 7.5mm

**C65-7280** Durham Tube, Neutral Glass, 25 x 6.5mm

**C65-7300** Durham Tube, Neutral Glass, 30 x 6.5mm

**C65-7320** Durham Tube, Neutral Glass, 35 x 7.5mm

**C65-7340** Durham Tube, Neutral Glass, 50 x 7.5mm

## TEST TUBES, POLYSTYRENE, NON-STERILE

A range of non-sterile, disposable polystyrene tubes.

**C65-8110** Round Tubes, 65 x 10mm. Pack of 6000

**C65-8130** Round Tubes, 75 x 12mm. Pack of 3500

**C65-8150** Round Tubes, 100 x 16mm. Pack of 1200

**C65-8170** Conical Tubes, 100 x 16mm. Pack of 1200

## TEST TUBES, POLYSTYRENE, STERILE

Sterile, disposable polystyrene tubes, with polyethylene screw cap.  
Pack of 450

**C65-8200** Test Tubes, Round End, 100 x 16mm

**C65-8220** Test Tubes, Conical End, 100 x 16mm



C65-6400 et seq



C65-6600 et seq



C65-8110 & C65-8130



C65-8220