

DISPLACEMENT VESSEL

Plastic body with metal base and copper overflow spout. Approx. size: 145 x 60mm.

P11-0500 Displacement Vessel

DISPLACEMENT VESSEL

Metal with overflow spout. Dimensions stated are approx. height x diameter.

P11-0610 Displacement Vessel 175 x 115mm

P11-0615 Displacement Vessel 100 x 50mm

PERSPEX CONTAINER

Approx. 50 x 50 x 90mm high. With line marked 80mm from base.

P11-0800 Perspex Container

PERSPEX CUBIC LITRE

100mm cube, holds one litre of water which weighs 1kg.

P11-0900 Perspex Cubic Litre

LITRE SET

The set consists of five different shaped clear plastic containers which each hold one litre of liquid.

P11-1000 Litre Set

GEOMETRICAL SHAPES SET, yellow plastic

Accurately made 3 dimensional shapes with size relationship of 50mm base section or diameter and 100mm vertical height.

P11-1201 Geometrical Shapes Set




P11-0500



P11-0610 & P11-0615



P11-0800



P11-0900



P11-1000



P11-1201