

CONDENSERS

CONDENSERS, BOROSILICATE GLASS

With push-on olives.

	Type	Effective length mm	Overall length mm
C30-1000	Liebig	250	440
C30-1050	Liebig	500	690
C30-1100	Allihn	250	440
C30-1150	Allihn	500	690
C30-1200	Graham	250	440



C30-1000

CONDUCTIVITY METERS

TOTAL DISSOLVED SOLIDS/CONDUCTIVITY METERS, STICK

Microprocessor based TDS and Conductivity meters with Automatic Temperature Compensation. Features large LCD, operation by a single button on the top, automatic switch off after 5 minutes of inactivity to conserve battery life.

Technical Data

	C30-4500	C30-4520	C30-4540
Range:	0-1999ppm	0-1999ppm 0.0-60°C	0-1999µS/cm
Resolution:	1ppm	1ppm; 0.1°C	1µS/cm
Accuracy:	±2% of full scale; 0.5°C		
ATC:	0-60°C (32-140°F)		
Environment:	0-50°C (32-122°F); 95% RH		
Dimensions:	180 x 50 x 25mm		
Weight:	80g		

Supplied complete with 2 x 1.5V batteries and full instructions.

C30-4500 Total Dissolved Solids Meter, 0-1999ppm

C30-4520 TDS/Temperature Meter, 0-1999ppm/0-60°C

C30-4540 Conductivity Meter, 0-1999µS/cm


Accessories

C30-4610 Calibration Solution, 1413µS/cm, 230ml

C30-4620 Calibration Solution Sachets, 1413µS/cm, 25 x 20ml

C30-4630 Calibration Solution, 1382ppm, 230ml

C30-4640 Calibration Solution Sachets, 1382ppm, 25 x 20ml



C30-4500