

FILTER CRUCIBLES, GLASS

All capacities and dimensions quoted are nominal

FILTER CRUCIBLES, PYREX GLASS

Gooch pattern with sintered disc.

Disc 20mm, capacity 10ml, Pyrex 3650 series. Pack of 5

	Porosity
C31-5000	1
C31-5050	2
C31-5100	3
C31-5150	4

Disc 30mm, capacity 30ml, Pyrex 3660 series. Pack of 5

	Porosity
C31-5350	1
C31-5400	2
C31-5450	3
C31-5500	4
C31-5550	5

Disc 40mm, capacity 60ml, Pyrex 3680 series. Pack of 5

	Porosity
C31-5750	1
C31-5800	2
C31-5850	3
C31-5900	4
C31-5950	5

ADAPTER FUNNEL, PYREX

For Pyrex filter crucibles. Stem o.d. 9mm, length 75mm. Pyrex 3690 series. Pack of 10. Rubber adapter **C31-6300** series required.

	For disc dia.	Open end i.d.	Height (exc. stem)
	mm	mm	mm
C31-6100	20	31	60
C31-6150	30	40	68
C31-6200	40	47	78

RUBBER ADAPTER

For filter crucibles.

	For disc dia., mm
C31-6300	20
C31-6350	30
C31-6400	40

FILTER CRUCIBLES, PORCELAIN

All capacities and dimensions quoted are nominal.

FILTER CRUCIBLES, PORCELAIN

Gooch pattern, with lid and perforated disc.

C31-6670 Capacity, 25ml



C31-5000 et seq



C31-6100 et seq and C31-6300 et seq in use



C31-6670