

MEMBRANE FILTERS

CELLULOSE NITRATE MEMBRANES, TYPE WCN, WHATMAN

For the majority of applications involving particles and cells. Suitable for all gravimetric and clarification techniques with aqueous solutions.

Plain white

	Pore size, μm	Dia., mm	Box qty
C40-2500	0.1	25	100
C40-2520	0.1	47	100
C40-2560	0.2	13	100
C40-2580	0.2	25	100
C40-2600	0.2	47	100
C40-2620	0.2	90	25
C40-2640	0.2	142	25
C40-2720	0.45	13	100
C40-2740	0.45	25	100
C40-2760	0.45	37	100
C40-2780	0.45	47	100
C40-2800	0.45	82	25
C40-2840	0.45	90	25
C40-2880	0.45	142	25
C40-2900	0.45	293	25
C40-2920	0.65	25	100
C40-2940	0.65	47	100
C40-2980	0.8	25	100
C40-3000	0.8	37	100
C40-3020	0.8	47	100
C40-3080	1.0	25	100
C40-3100	1.0	47	100
C40-3140	3.0	25	100
C40-3160	3.0	47	100
C40-3200	5.0	25	100
C40-3220	5.0	47	100
C40-3240	5.0	90	25

MEMBRANE FILTER HOLDERS

FILTER HOLDERS, VACUUM TYPE, WHATMAN

With borosilicate glass funnel and base, metal clamp and choice of sintered glass or stainless steel membrane support. Complete with gasket, 'O' ring and silicone rubber stopper on outlet for fitting into filter flask. 25mm holders have No. 5 stoppers fitting 125ml flasks. 47mm holders have No. 8 stoppers fitting standard manifolds and 1 litre flasks.

25mm diameter

	Support	Reservoir, ml	Prefilter, mm
C40-4200	Sintered glass	25	16
C40-4250	Sintered glass	50	Not required
C40-4300	Stainless steel	25	16

47mm diameter

	Support	Reservoir, ml	Prefilter, mm
C40-4350	Sintered glass	300	35
C40-4400	Stainless steel	300	35

FILTRATION MANIFOLDS, WHATMAN

For use with vacuum type filter holders. Stainless steel with individual PTFE stopcocks.

C40-4600 3 Port Model



C40-4200 et seq,
shown in use
with C40-4600