

FILTER HOLDERS, SYRINGE TYPE

Polycarbonate, autoclavable, luer lock end fittings.

C40-4700 Filter Holder, for filter diameter 13mm. Pack of 12

C40-4720 Filter Holder, for filter diameter 25mm. Pack of 12

FILTER HOLDERS, SYRINGE TYPE, WHATMAN

Stainless steel, with Luer lock fittings to connect to a standard syringe. PTFE gaskets and 'O' ring permit autoclaving of the membrane in place without the filter sticking to the holder.

C40-4800 Filter Holder, for filter diameter 13mm

C40-4850 Filter Holder, for filter diameter 25mm

IN-LINE FILTER HOLDER

Polycarbonate, autoclavable, tubing connectors on each end, for use with 9 - 13mm I.D. tubing. Holds 47 or 50mm filter membranes.

C40-4900 In-Line Filter Holder, 47mm. Pack of 5

FILTER HOLDER WITH RECEIVER, NALGENE

Polysulfone, autoclavable. Features graduated receiver, two side-arms enable connection to vacuum line, tubing adapter to accept 6 to 8mm I.D. tubing. Supplied with receiver cover to allow storage of the sterile filtrate. Holds 47 or 50mm filter membranes.

	Upper chamber capacity, ml	Receiver capacity,ml
C40-4960	250	250
C40-4970	500	500
C40-4980	500	1000

DISPOSABLE FILTERS

SYRINGE FILTERS, SARTORIUS, MINISART

Suitable for small volumes. Comprises a cellulose acetate membrane sealed into a clear plastic housing with Luer lock fittings on inlet and outlet. Filter area 5.3cm². Will filter in either direction. Supplied sterile, individually sealed.

C40-5250 Syringe Filter, 0.2µm pore size. Pack of 50

C40-5300 Syringe Filter, 0.45µm pore size. Pack of 50

STERILISING FILTER UNITS, NALGENE

Pre-sterilised disposable single use filter units, polystyrene housing, cellulose nitrate membrane, pore size 0.2µm. Supplied sterile, individually packed, outer pack of 12.

	Upper chamber capacity, ml	Receiver capacity,ml
C40-5600	150	150
C40-5610	250	250
C40-5620	500	500
C40-5630	500	1000

Filter Funnels, see Funnels section

Filter Pumps, see Pumps section

First Aid Kits, see Safety section



C40-4960 et seq



C40-5250



C40-5600