

**ABSORPTION****ABSORPTION TOWER, Glass**

With ground neck and ground inlet. Without stopcock and tubing adapter. Dimensions stated are overall height x diameter. Capacity, 250ml.

**C10-0200** 315 x 45mm, pack of 8

**ABSORPTION TUBES, Glass**

Straight form, with one bulb and tube for bending. Dimensions stated are overall length, excluding narrow tube, x diameter.

**C10-0400** 100 x 12mm

**C10-0450** 150 x 20mm

**C10-0500** 200 x 25mm

**ABSORPTION TUBES, Glass**

Plain U- form. Dimensions stated are overall length x diameter.

**C10-0550** 100 x 12mm

**C10-0600** 125 x 15mm

**C10-0650** 150 x 20mm

**ABSORPTION TUBES, Glass**

U-form, with side arms. Dimensions stated are overall length x diameter.

**C10-0700** 100 x 12mm

**C10-0750** 150 x 20mm

**ABSORPTION TUBES, Glass**

U-form, with side arms and ground-in stopper-stopcocks. Dimensions stated are overall length x diameter.

**C10-0800** 100 x 12mm

**C10-0850** 150 x 20mm



C10-0200



C10-0550 et seq



C10-0700 et seq



C10-0800 et seq