

ADAPTERS

ADAPTERS, Glass

T-shape. Dimensions stated are external diameter of limb.

C10-1000 5mm, pack of 5

C10-1050 7mm, pack of 5

ADAPTER, Glass

T-shape with two stopcocks, with limbs 7mm external diameter.

C10-1100 Adapter with stopcocks

ADAPTERS, Polypropylene

Straight form, with tapered stepped ends. Dimensions stated are min/max respectively for each end of connector.

C10-1200 4/7 & 4/7mm o.d., pack of 10

C10-1250 4/7 & 7/12mm o.d., pack of 10

C10-1300 7/12 & 7/12mm o.d., pack of 10

C10-1350 4/7 & 10/15mm o.d., pack of 10

C10-1400 7/12 & 10/15mm o.d., pack of 10

ADAPTERS, Polypropylene

Straight form, two part push-fit connectors. Dimensions stated are bores of tubing accepted.

C10-1510 6 to 8mm, pack of 10

C10-1560 8 to 10mm, pack of 10

C10-1600 10 to 12mm

C10-1650 12 to 14mm

ADAPTERS, Polypropylene

T-shape, with tubing grip serrations, autoclavable. Dimensions stated are bores of tubing accepted.

C10-1810 6mm, pack of 10

C10-1860 8mm, pack of 10

C10-1910 10mm, pack of 10

C10-1960 12mm, pack of 10

C10-2010 14mm, pack of 10

ADAPTERS, Polypropylene

Y-shape, with tubing grip serrations, autoclavable. Dimensions stated are bores of tubing accepted.

C10-2110 8mm, pack of 10

C10-2160 10mm, pack of 10

C10-2210 12mm, pack of 10

C10-2260 14mm, pack of 10

ADAPTER VALVES, Polypropylene

Straight form, two part push-fit connector incorporating a non-return valve. Dimensions stated are bores of tubing accepted.

C10-2410 8 to 10mm, pack of 10

C10-2470 12 to 15mm, pack of 10




C10-1000 et seq



C10-1200 et seq



C10-1510 et seq



C10-1810 et seq



C10-2110 et seq



C10-2410 et seq