

AEROSOLS

ADHESIVE SPRAY

Produces a fine mist of adhesive, suitable for many substrates. Particularly useful for paper, card and film, temporary positioning of small components, allows short term repositioning. Not suitable for styrene based products.

C10-3010 Adhesive Spray, 400ml

ANTISTATIC SPRAY

Prevents the build-up of static, thus eliminating dust and unpleasant shocks and reducing fire hazards in areas such as computer rooms, etc. Supplied in a pump action bottle.

C10-3050 Antistatic Spray, 400ml

FREEZER SPRAY

For rapid cooling of electrode components to locate faults or to prevent heat damage during soldering. Also useful for shrink fitting metal parts and restoring broken mercury threads in thermometers. Reduces temperature down to -50°C. Safe on most materials.

C10-3060 Freezer Spray, 200ml

LUBRICANT SPRAY, WD-40

General lubricant, penetrator, anti-moisture agent and cleaning oil.

C10-3080 Lubricant Spray, WD-40, 200ml

PTFE SPRAY

Coats equipment with a layer of chemically inert, anti-stick PTFE to provide self lubricating surfaces. Suitable for glass joints, bearings, gaskets, valve parts, stopcocks, gears etc. It is also an excellent mould release agent.

C10-3110 PTFE Spray, 200ml

SILICONE GREASE SPRAY

Gives a thin but uniform coating of silicone grease, is oxidation resistant and water repellent. Uses include coating interchangeable glass joints, waterproofing electrical circuits, lubricating moving parts and preserving rubber.

C10-3120 Silicone Grease Spray, 400ml



C10-3010



C10-3050



C10-3060



C10-3080



C10-3110



C10-3120