

ALUMINIUM BLOCKS

For use with all Stuart block heaters. Made from anodised aluminium. Each block features an extra hole for a thermometer if required. All blocks (except **C14-6518** and **C14-6519**) measure 75 x 95 x 50mm (wxdxh) and can be used in any combination. **C14-6518** and **C14-6519** measure 150 x 95 x 61mm (wxdxh).

	Tube (dia./type)	No. of holes	Hole size, mm (dia. x depth)
C14-6501	Plain block	-	-
C14-6502	10mm	20	10.5 x 47
C14-6503	12mm	20	12.6 x 47
C14-6504	12mm	20	12.6 x 33
C14-6505	13mm	20	13.5 x 47
C14-6506	16mm	12	16.5 x 47
C14-6507	19mm	8	19.5 x 47
C14-6508	2ml tubes	20	10.5 x 33
C14-6510	1.5ml tubes	20	10.7 x 14
C14-6511	25mm	6	25.5 x 47
C14-6512	28mm	6	28.0 x 47
C14-6513	0.5ml tubes	30	7.8
C14-6514	30mm	4	30.1 x 47
C14-6515	33mm	4	33.2 x 47
C14-6516	0.2ml tubes	48	6.1
C14-6517	0.2ml strip	10 x 8	-
C14-6518	96 well plate	-	-
C14-6519	384 well plate	-	-

SAFETY COVERS

Tough Perspex cover, easily fitted without tools. Protects the user from 'spitting' tubes.

C14-6550 Safety cover for 2 block models

C14-6551 Safety cover for 3 block and dual control models



C14-6501 et seq



C14-6550